IMPIANTI DENTALI "IS-TWO IS-FOUR IS-FIVE CYROTH KAROTH SLC TAC & TAC Zygo"

DESCRIZIONE

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DENTAL IMPLANTS "IS-TWO IS-FOUR IS-FIVE CYROTH KAROTH SLC TAC & TAC Zygo"

DESCRIPTION

AoN Implants dental implants are endosseous screws designed and sized for all applications in the field of odontological surgery. They can be applied in previously prepared sites in edentulous areas in accordance with an operating sequence described in apposite surgical techniques

INTENDED USE

The dental implants are designed for applications in odontological implantology.
The characteristic indications of use are:
single edentulous
partial or total edentulous
distal edentulous
post-extraction replacement
All the implants are mono-use and must not be reused after eventual removal.

MATERIALS

All materials are of high quality and conform to reference standards. AON dental implants are realised in titanium alloy, supplied by qualified companies and conforming to international reference standards for titanium that can be used in implantology.

The titanium Grade 4 guarantee an high degree of biocompatibility

The dental implants, after cleaning and decontamination, are sold sterilised with gamma rays inside of the final packaging and ready for its intended use.

The dental implants are sold in single units in small double plastic bottles designed with such shape and transparency as to immediately identify the contents which, in turn, are placed inside a secondary packaging box with external labelling indicating its contents. The sterility is guaranteed up until breakage or damage of the outer plastic wrapping. The double packaging, therefore, assures the maintenance of sterile conditions until the removal for the placement in the oral cavity. It is recommended that the opening of the sterile package always be done in an aseptic environment and possibly over worktops with disinfected or sterile surfaces.

CONTRAINDICATIONS

The use of implant therapy is clearly not recommended for the patients in unfit conditions for a general surgery. The device is not indicated for use in patients with reduced immunity, coagulation problems, endocrine disorders, severe rheumatic or bone disease, liver cirrhosis or other serious illnesses. Other contraindications are related to:

- Patients with aesthetic expectations not realizable
- Psychotic patients and/or with emotional instability
- Patients submitted to local radiotherapy
- Patients with severe bruxism and/or serious cranio-facial and maxilla-mandibular discrepancies
- insufficient quality of bone tissue and/or quantity of the bone
- systemic diseases of the blood
- forms of diabetes - not controlled
- cardiovascular problems
- infarct pathology
- metabolic disorders
- immunodepression
- bone growth - not mature
- ongoing inflammatory process
- bone tumours
- bone tissue loss that would prevent primary stabilisation
- maxillary-mandibular neoplasms
- radiation therapy in the maxillo-facial area
- mucosal pathology in general
- Smokers of more than 10 cigarettes per day

Moreover, the use of implants is not recommended in patients who take drugs, abuse of alcohol and smoke, or patients with poor oral hygiene, inasmuch such behaviour in general can compromise the success of surgical treatment.

USE BY DATE

There are no special storage conditions although it is always preferable to store in a cool, and dry place.
The expiration period of 5 years is from the date of sterilisation which is always indicated on the label applied to the external packaging and in the bottle within the package. Disposal of waste after use must follow the existing procedures for the disposal of medical waste products.

INFORMATION FOR THE PATIENT

The law provides that the patient be properly informed by the doctor concerning the risks and surgical complications that could arise. The doctor must, moreover, instruct the patient as he always maintain the visits of control and return to his surgery whenever unforeseen situations relating to the implant occur.

Use in an outpatient setting is reserved for qualified and authorized professionals in accordance with the law, specialists in oral surgery and prosthetic reconstruction. The devices are used on adult patients who have completed their craniofacial development, with no age restrictions.

TRACKING

The AON Dental implants have labels which allow traceability of the product, pursuant to current regulations. These labels will be pasted: one on the patient's medical records archived in the surgery, the other is pasted in the "Implant Identity Card" document that will be delivered to the patient.

SIDE EFFECTS

These devices should not be used in allergic patients or sensitive to titanium. Other side effects related to the application of the implants are:

- hematoma of the effected region
- acute infections and late infections in the region operated on
- temporary or persistent functional modification of the nerves in the affected anatomical region
- vein thrombosis, pulmonary embolism, cardiac arrest
- modification of position and/or mobilisation of the implant
- reduction of bone tissue due to crestal reabsorption
- relaxation of the gingival tissue
- pathological bone fractures caused by modifications in load
- allergic or metallosis reactions in the regions adjacent to the implant

INSTRUCTIONS AND MODALITY OF USE

The utilisation of IS-one dental implants is exclusively reserved to highly qualified odontological medical personnel with appropriate training and experience in implantsurgery. Before the application the following tests are suggested: Interview with the patient in order to verify the needs and the expectations.

- Clinical examination
- Rx panoramic
- Rx dental scan or Rx 3D cone beam
- If the traditional radiology is not enough to plan the operation it is appropriate to design it on the basis of the TAC results with dedicated simulation program
- Study models and diagnostic waxing
- Surgery dima
- 1) Preparation of the implant site according to the sequence envisioned in the surgical technique
- 2) Select the appropriate screw, as to type and size
- 3) Extract the small bottles from the box.
- 4) Remove the seal and extract the bottle inside
- 5) Unscrew the cap thereby providing access to the screw (there are two variants: implant with connection without hexagon or implant with pre-assembled mounter)
- 6) Take the screw from the bottle with the apposite instrument and position it in the opportunely prepared site
- 7) Proceed with the insertion by threading it in, without exerting excessive pressure
- 8) with threaded cap (placed in the lower part of the bottle) and suture the gingival flaps .
- The pre-assembled mounter can also serve as a transfer or temporary abutment.
- 9) Wait the time needed for bone integration
- 10) On completion of osseointegration, complete the incision of the gingival tissue located above to expose the head of the implant.
- 11) Remove the closure screw
- 12) Apply the healing screw and suture the soft tissues around the prosthetic component.
- In any case, it is always compulsory the knowledge of the specific protocol that must be used reported to the product in use.

WARNINGS

An early loss of anchorage might be verified after insertion of implants. The potential causes are:

- lack of quantity and/or quality of bone
- Poor oral hygiene and/or genetic diseases
- allergic and/or hypersensitivity reactions to the metal
- Disturbances due to the presence of the device
- inadequate healing
- Peri-implantitis due to gingival irritation

The risks of failure are related to an use which is not adequate to the anticipated prescriptions and also to an incorrect evaluation of the patient prior to application. In fact: visits, assessments of the quality and quantity of bone are necessary. Panoramic radiograph performed, localised, and other types of examinations, not just radiographic can be necessary.

ACCIDENT REPORTS

Users and/or patients must report to the manufacturer AON IMPLANTS and to the competent authority of the member state in which they are established, any serious accident that has occurred in relation to the device.

The Manufacturer AoN Implants has drawn up the Summary document on Safety and Clinical Performance (SSCP) which is available on request and made public through the EUDAMED platform.

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