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LEGENDA SIMBOLI

SURGICAL INSTRUMENTS AND ACCESSORIES

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DESCRIPTION

Instruments and Accessories include: Contra-angle fixing keys, Drills, Final Drills, Sterile Final Drills, Guides, Lance Drills, Tapping Devices, Mucotomes, Countersinks, Extensions, Probes, and Stops.

Among the accessories, zygomatic demo implants are also available for the correct selection of implant size by inserting them into the prepared site.

Instruments can be delivered pre-assembled in sets in specific sterilisation trays.

INTENDED USE

The Instruments and Accessories are specific devices for oral implantology, used in the preparation of the implant site for the subsequent placement of a dental implant.

They must be fit on equipment designed for this purpose (handpiece – micromotor) and can only be used by qualified personnel (dentists or surgeons).

With regard to the limits of use, please bear in mind that all items must be sterilised before use.

A line of “final” single-use drills specifically designed for Is-Five implants is available, marketed in sterile packaging and supplied in combination with this implant.

All AON Implants products are identified with a UDI carrier on the packaging, while instruments are marked with the batch number.

MATERIALS

Instruments and Accessories are made of AISI 420 and AISI 630 stainless steel (the latter also known as 17-4 PH stainless steel).

The devices undergo heat treatment hardening and are always electropolished with a specific acid bath followed by validated ultrasonic cleaning.

The devices have physical and mechanical properties suitable for use, according to the regulations in force for the manufacturing of surgical instruments.

The zygomatic demo implant accessories are made of biocompatible grade 4 titanium.

CONTRAINDICATIONS

The instruments have been designed for the preparation of the implant site. As regards the contraindications of implant therapy, please refer to the specific instructions for use provided.

Moreover, the use of implants is not recommended in patients who take drugs, abuse of alcohol and smoke, or patients with poor oral hygiene, inasmuch such behaviour in general can compromise the success of surgical treatment.

STORAGE, MAINTENANCE AND DISPOSAL

Store the instruments in a dry place and do not expose them to heat sources or sunlight. The instruments must be protected to heat sources or sunlight. The instruments must be protected from corrosive agents.

Before each use, check the degree of wear of drills and that they do not show mechanical damage; if in doubt, do not use them. Before disposal by conventional means provided for medical waste, used drills must be cleaned and sterilized.

INFORMATION FOR THE PATIENT

The law provides that the patient be properly informed by the doctor concerning the risks and surgical complications that could arise.

PATIENT SELECTION CRITERIA

The manufacturer AON IMPLANTS does not place any restrictions in terms of patient selection for the use of Instruments and Accessories, as they are made with biocompatible materials suitable for their intended use and can be used on adult patients who have completed craniofacial growth, without age limits. However, the dental specialist is solely responsible for the use and selection of such devices.

SPECIFIC RESPONSIBILITIES

The AON Implants Instruments and Accessories for oral implantology meet the general safety and performance requirements set out by Regulation (EU) 2017/745-MDR. Any use of the products other than that specified is to be considered as “improper use”, releasing the manufacturer from any liability.

WARNINGS

An early loss of anchorage might be verified after insertion of implants. The potential causes are:

- lack of quantity and/or quality of bone
- Poor oral hygiene and/or genetic diseases
- allergic and/or hypersensitivity reactions to the metal
- Disturbances due to the presence of the device
- inadequate healing
- Peri-implantitis due to gingival irritation

The risks of failure are related to an use which is not adequate to the anticipated prescriptions and also to an incorrect evaluation of the patient prior to application. In fact: visits, assessments of the quality and quantity of bone are necessary. Panoramic radiograph performed, localised, and other types of examinations, not just radiographic can be necessary.

TRACKING

The AON Dental implants have labels which allow traceability of the product, pursuant to current regulations. These labels will be pasted: one on the patient’s medical records archived in the surgery, the other is pasted in the “Implant Identity Card” document that will be delivered to the patient.

INCIDENT REPORTS

Users and/or patients must report to the manufacturer AON IMPLANTS and to the competent authority of the member state in which they are established, any serious accident that has occurred in relation to the device.

The Manufacturer AoN Implants has drawn up the Summary document on Safety and Clinical Performance (SSCP) which is available on request and made public through the EUDAMED platform.

ADDITIONAL INFORMATION

The Instruments and Accessories must always be sterilised before use.

The devices are supplied in NON-STERILE condition. Only the Is-Five final drill is marketed in single-use sterile packaging.

The Instruments and Accessories are NOT devices intended for use by non-experts.

SIDE EFFECTS

The Zygo Demo device is made of titanium and therefore should not be used for patients who are allergic or sensitive to this material.

For side effects of implant therapy, please refer to the specific instructions for use provided.

METHOD OF USE

AON Implants S.r.l. Instruments and Accessories must be used exclusively for the use indicated in the product catalogue. Different use could cause damage and the consequent loss of the original features and functionalities or breakage with consequent risks of injury.

Method of use for rotating instruments:

Ensure that rotating instruments are perfectly secured to the micromotor before starting it. Also, ensure not to exceed the maximum rpm and the number of uses indicated in the limitations of use.

The multi-use final drills, single-use drills, drills, and extensions must never be used at more than 800 revolutions per minute; the only exception is lance drills, which must not exceed 1,000 revolutions per minute and must never be used without a saline cooling solution to prevent bone overheating and possible necrosis that may ensue.

If a BHC protocol without irrigation is used, as indicated in the product catalogue and surgical protocols, all drills must never exceed 70 revolutions per minute.

Lance drills cannot be used in a BHC protocol and must therefore be used up to a maximum of 1,000 revolutions per minute always with the aid of irrigation with cold saline solution for proper cooling.

If necessary, the user must plan breaks during surgery to allow the bone and the drill to cool naturally.

The depth to be reached is indicated by the reference notches on the drills or coincides with the end of the cutting edge.

The sequence of drills to be used for each type of implant and any specific indications are reported in the surgical protocol.

It is necessary to check the sharpness of the instrument and not exceed the number of uses indicated in the product description. Multi-use drills can be used up to a maximum of 30 times.

In case of doubt about the sharpness, the drill must always be replaced.

The Is-Five final drills are single-use and must not be reused to avoid cross-contamination and loss of performance.

The contra-angle mounter keys for implant placement and the prosthetic keys for tightening the prosthetic screws must be used at a maximum of 10 revolutions per minute.

The maximum recommended torque for prosthetic keys is 40 Ncm.

The maximum recommended torque for mounter keys is 70 Ncm.

CLEANING AND STERILISATION

The cleaning and sterilisation processes aim to safeguard the health of patients and all people working at the clinic. If the cleaning and disinfection processes are carried out correctly, they will enhance the quality of the tools and prolong their efficiency over time.

It is advised not to perform unnecessary sterilisation cycles on drills.

Before sterilisation, drills must always be thoroughly cleaned. Cleaning

Cleaning can be carried out manually with warm water and a suitable non-corrosive neutral detergent, using plastic or nylon brushes (never steel wool or metal brushes) to remove any organic residues. Soft brushes do not damage the cutting edges but allow for the complete removal of residues.

Therefore, manual washing is effective only if carried out with utmost care. Be careful of the risk of injury, especially with drills not yet sterilised!

Avoid using excessively aggressive detergents, such as hydrogen peroxide or high concentrations of sodium hypochlorite, or other solutions containing high concentrations of chlorine.

In any case, when using the chosen detergent, it is necessary to strictly follow the manufacturer’s instructions. Cleaning in ultrasonic baths is certainly safer for the operator compared to manual cleaning. It is effective and saves time.

When using ultrasonic baths, it is always necessary to check that the instruments are isolated during treatment, to avoid contact between pieces that may damage the cutting surfaces.

It is always necessary to individually check the devices to ensure proper cleaning. Suitable non-corrosive neutral detergents (preferably containing a corrosion inhibitor) can be used, always in accordance with the manufacturer’s instructions for the detergent and ultrasonic device used.

Dip the drills completely and leave them in the washing liquid for the time recommended by the ultrasonic device manufacturer. Change the washing solution at the intervals recommended by the manufacturer. After ultrasonic washing, thoroughly rinse the drills to remove the washing solution.

Immediately after rinsing, before sterilisation, always inspect and dry them.

Do not the devices wet; after possible rinsing with “hard” water, dipping in ethyl alcohol for additional sanitising and dehydration can be carried out.

Always provide for drying even before sterilisation when using the relevant sterilisation trays.

If autoclaved in a suitable bag for steam sterilisation, it is recommended not to package products made of different material in the same bag, in order to prevent galvanic corrosion phenomena.

Sterilisation

First of all, remember that sterilisation does not replace cleaning.

Sterilisation is a process that kills all microorganisms and must be carried out with utmost care after cleaning. The most common and effective sterilisation methods that may be verified at a clinic with “Spore Testing” are as follows:

Autoclave: the tools are treated with high-pressure steam. Recommended standard time*: 20 minutes at a temperature of 120°C (250°F). Sterilising liquid: distilled water.

Chemical vapours: the tools are treated with high-pressure chemical vapours.

It is strictly necessary to dry the tools before starting sterilisation. Recommended standard time *: 20 minutes at a temperature of 130°C (270°F).

Sterilising liquid: special solutions.

Dry heat sterilisation: sterilisation of the tools takes place by hot air. It is strictly necessary to dry the tools before starting sterilisation.

Recommended standard time: 60 to 90 minutes at a temperature of 150°/170°C (300°/320°F).

The recommended time and temperatures may vary depending on the nature and load of your equipment and do not include warm-up time.

Always follow, however, the instructions provided by the equipment’s manufacturer.

For all reusable instruments, AON IMPLANTS recommends sterilisation in an autoclave in accordance with a validated method in compliance with standard UNI EN ISO 17665-1.

Only the Is-Five final drill is marketed already sterile in single packaging and is single-use.

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